



**Vaccine Storage Unit Temperature Log** Month/Year \_\_\_\_\_ Storage Unit \_\_\_\_\_ Provider \_\_\_\_\_ VFC# \_\_\_\_\_


- 1) Check storage unit temperature twice a day. Record the temperature by putting an X in the corresponding box and filling in all required information in the column for that day.
- 2) If using a Data Logger, check that the LED light is green. Put an X in the "LED Green" box for the corresponding day and time.
- 3) If temperatures are in the shaded zone or the Data Logger LED light is red, DO NOT USE or DISCARD the vaccine. Obtain a Vaccine Incident Report ([www.immunization.mt.gov](http://www.immunization.mt.gov) under "VFC" link) and follow the instructions on the top of page 1. Call the Montana Immunization Program at 444-5580. **VFC temperature logs must be saved for three years.**

| Day of Month            | 16   |   | 17 |    | 18 |    | 19 |    | 20 |    | 21 |    | 22 |    | 23 |    | 24 |    | 25 |    | 26 |    | 27 |    | 28 |    | 29 |    | 30 |    | 31 |  |
|-------------------------|------|---|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|--|
| Time of Day             | am   | pm  | am | pm | am | pm | am | pm | am | pm | am | pm | am | pm | am | pm | am | pm | am | pm | am | pm | am | pm | am | pm | am | pm | am | pm |    |  |
| Staff Initials          |      |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |  |
| LED Green               |      |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |  |
| <b>Refrigerator (F)</b> | ≥49° | <b>DO NOT USE vaccine that has been stored at these temperatures (shaded grey). See instructions (Step 3) at top of page.</b> |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |  |
|                         | 48°  |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |  |
|                         | 47°  |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |  |
|                         | 46°  |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |  |
|                         | 45°  |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |  |
|                         | 44°  |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |  |
|                         | 43°  |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |  |
|                         | 42°  |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |  |
|                         | 41°  |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |  |
|                         | 40°  |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |  |
|                         | 39°  |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |  |
|                         | 38°  |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |  |
|                         | 37°  |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |  |
|                         | 36°  |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |  |
|                         | 35°  |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |  |
|                         | 34°  | <b>DO NOT USE vaccine that has been stored at these temperatures (shaded grey). See instructions (Step 3) at top of page.</b> |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |  |
| 33°                     |      |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |  |
| 32°                     |      |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |  |

|                    |     |   |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
|--------------------|-----|---|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|
| LED Green          |     |   |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| <b>Freezer (F)</b> | >8° | <b>DO NOT USE vaccine that has been stored at these temperatures (shaded grey). See instructions (Step 3) at top of page.</b> |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
|                    | 7°  |   |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
|                    | 6°  |   |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
|                    | 5°  |   |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
|                    | 4°  |   |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| 3°                 |     |   |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |

# Vaccine Storage Unit Trouble-Shooting Log

Record issues with storage unit and corrective actions taken in the table below.  
 (e.g., near out-of-range temperatures, malfunctions, repairs)



| Date | Time | Description of Issue | Corrective Action Taken | Staff Initials |
|------|------|----------------------|-------------------------|----------------|
|      |      |                      |                         |                |
|      |      |                      |                         |                |
|      |      |                      |                         |                |
|      |      |                      |                         |                |
|      |      |                      |                         |                |
|      |      |                      |                         |                |
|      |      |                      |                         |                |